

Airwave Therapy
Evolution in Zero Pressure

ROBER

Distributed by:

www.roberlimited.com

Tel. +44 (0) 1246 245404
@ info@roberlimited.com
@RoberLtd



CE
EN60601-1



NoDec®

© 2015 NoDec® is a registered trademarks of Rober Limited

Rober Limited, Wheatbridge Road, Chesterfield, Derbyshire, S40 2AB, UK.

ROBER

Rober Ltd reserves the right to make technical changes.

Zero Pressure Mattresses

NoDec®
PRODUCT SELECTOR

Patient Mobility	 Limited	 Limited	 Limited	 Immobile	 Intensive Care	 Can't be turned
Stage	 Stage 1-4	 Stage 1-4	 Stage 1-4 (Heal Ulcers)	 Stage 1-4	 Stage 1-4	 Stage 1-4
Risk Category	 High	 High - Very High	 High - Very High	 Very High	 Very high	 High/Very high
Recommended Repositioning Frequency*	 2-4 hrs	 4-6 hrs	 4-6 hrs	 8-12 hrs	 8-12 hrs	 6-8 hrs
Cell Cycle	 2 Cell	 2 Cell	 2 Cell	 3 cell	 3 cell	 3 cell
Weight	 30-180 kg	 20-220 kg	 20-220 kg	 20-250 kg	 20-250 kg	 30-400 kg
The best NoDec® Product for you	 AIRSTREAM Overlay Mattress <small>(Other Mattress variations are available)</small>	 AIRSTREAM2 Mattress	 AIRSTREAM THP (Total Heel Protection) Mattress	 AIRSTREAM3 Mattress <small>(Other Mattress variations are available)</small>	 WIZARD Mattress	 BARI Mattress

*International Guidelines (NPUAP/EPUAP/Pan Pacific) states 'repositioning frequency should be influenced by the support surface used. (Strength of Evidence = A). The repositioning frequency should depend on the pressure-redistributing qualities of the support surface.